


## QC Benchmark: a streamlined web application to comprehensively evaluate instrument performance and direct troubleshooting

Benjamin A. Neely  
Paul E. Anderson  
W. Clay Davis  
Magnus Palmblad

*Identification of certain commercial equipment, instruments, software or materials does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the products identified are necessarily the best available for the purpose.*



# System Suitability – how you test is important





## Metrics

Track speed  
High-speed corners

Track speed  
Low-speed corners  
High-speed corners  
Braking  
etc.




## Modern Mass Spectrometer




## Metrics

22k MS/MS  
3552 PSMs  
103 Peptide groups  
77% coverage of BSA



75k MS/MS  
28,469 PSMs  
21,689 Peptide groups  
4140 Protein groups

# Proteomics - System suitability sample and metrics



## Complex Digest:

HeLa  
Yeast  
*E. coli*  
*C. elegans*  
Arabidopsis

## Basic metrics:

Ex.  
75k MS/MS  
28,469 PSMs  
21,689 Peptide groups  
4140 Protein groups

## ID-based Metrics:

- Detection of specific peptides
- Rates of peptide identification
- Coverage of known proteins


## ID-free Metrics:

- Peak width/shape
- TIC
- Mass error
- Charge state distribution
- ...
- etc.

- Determine whether system is working  
... from LC to MS/MS
- Communicate quality across time and labs

## QC – the current toolset

### ID-Free + ID-based



### NIST MSQC Pipeline



**Attention:** Support for NIST MSQC Pipeline is discontinued

Tools derived/inspired by MSQC:  
MassQC  
Metriculator  
QuaMeter

### Specific application QC tools

- TMT - TVT: Triple Knockout Proteome Standard Visualization Tool

### ID-Free only

- spRoCoP - Statistical Process Control in Proteomics
- Raw Beans
- QC-ART


- Local and/or web-based versions
- Can work on any lab's QC sample
- Some require the presence of predefined  $m/z$
- Can be used for quality controls within experimental runs

*\* this is not an exhaustive list, please excuse any omissions.*

# Yet another tool!!!

Focus on a weekly QC sample run (known sample, known conditions) for system suitability, so that metrics/reports can be made **without user specified parameters**.

Primary User




TIC  
Peak Widths  
PSMs




Universal material and  
methods makes quality  
comparable and more useful

Core Lab




Report pre-run  
performance to users  
Report instrument  
performance to  
future users

# QC Benchmarker – Goals and Rationale




- 1) Use commercial digest standards and standard methods
- 2) Upload raw file to web-app
- 3) Return high-level summary (good, marginal, bad) with in-depth reports for troubleshooting

# QC Benchmark Implementation – Focus on automation and speed



- Website:
- Containerized Linux OS with R, Python, Redis and Celery
  - TPP (wine) used for conversion and identifications
  - Existing R packages: mzR, MSnbase, pepXMLTab

ID-free + ID-based metrics




	Sample	LC	Source	MS1	MS2
ID-free	amount loaded on column	Peak-widths, peak capacity	stability	TIC	Total signal
ID-based	% spectra matched	pH of mobile phases	z-dist of PSMs	mass error	mass error; # PSMs

*\*there are more metrics than these*





# Development of QC Benchmark




The development of data used for QC Benchmark testing required using a standard sample, gradient and MS/MS. For now this has focused on thermo instruments and a 2 hour 200ng HeLa injection.

- Generated ~70 runs specifically looking at:
- Mobile phase pH (0.1% formic acid to no formic acid)
  - Different amounts of HeLa on column (50 to 1000 ng)
  - Sub-optimal collision energy
  - Attempted to degrade HeLa
  - Runs before and after MS cleaning
  - Attempted to decrease transfer efficiency
  - Mass error (couldn't get over 8 ppm)
  - Running on a different platform (QE classic v Lumos)
  - 15 runs across a year



# Example of a performance threshold – PSMs



- Thresholds are based on instrument and type of acquisition
- Since each instrument is unique, qualitative thresholds may not be applicable, but metrics are still valid






# Machine Learning-based prediction of mobile phase pH and pump performance





- Decreased correlation may indicate mobile phase pH is off
- A deviation of the intercept indicates the pumps aren't performing correctly


# Autocorrelation-based peak width determination



# Autocorrelation-based peak width determination






30 second peak at base  
FWHM = 11.7 sec



FWHM




# Autocorrelation-based peak width determination

- Perform across the whole gradient






Mean FWHM

# Autocorrelation-based measure of spray stability



Shifted all points by 50 ms



# Autocorrelation-based measure of spray stability



Shifted all points by 50 ms




## Autocorrelation-based measure of spray stability



- The lag-1 autocorrelation generates a summary of peak stability across the run
- Still working to determine what constitutes good, marginal or bad

ID-free + ID-based metrics



	Sample	LC	Source	MS1	MS2
ID-free	amount loaded on column	peak widths, peak capacity	stability	TIC	Total signal
ID-based	% spectra matched	pH of mobile phases	z-dist of PSMs	mass error	mass error; # PSMs

*\*there are more metrics than these*

---

msqc.live

Documentation

About

Disclaimer

Privacy Concerns

## QC Benchmark **BETA**

Save this url for future reference: [http://34.74.220.180/index/qc\\_benchmarker\\_lbpqsk56](http://34.74.220.180/index/qc_benchmarker_lbpqsk56)

Drop raw file here

By uploading a file you acknowledge reading the disclaimer

---

msqc.live

Documentation

About

Disclaimer


Privacy Concerns

## QC Benchmarker **BETA**


Save this url for future reference: [http://104.196.70.201/index/qc\\_benchmarker\\_lbpqsk56](http://104.196.70.201/index/qc_benchmarker_lbpqsk56)

---


2018-10-18\_HeLa\_OTOT



2018\_06\_27\_HeLa\_OTOT



2018-05-04\_HeLa\_OTOT



*\* click on colored  
modules for full reports  
underneath*

## Ongoing development (v1)

- Summary statistics and figures related to more “classic” quality measures
- Test on non-Thermo platforms
- Generate data and test for Water’s *E. coli* and Promega Yeast
  - please recommend any commercially available digest
- Establish relevant thresholds per platform
- Allow for shorter gradients or different loading amounts
- Model column heater effects and more pH data points

## Future development (v2)

- Allow user to specify thresholds
- Allow user to import specific reference file (ex. experiment specific conditions and QC)
- Locally evaluate temporal variations on large data sets
- Evaluate mass spec-based metabolomics data

---

# Questions

---